

Evaluation of run controls for cobas 6800 MPX and HEV assay

Marco Koppelman

Head of Nucleic Acid testing (NAT) laboratory,
National Screening laboratory of Sanquin (NSS) in the
Netherlands

Satellite Meeting before IPFA-PEI 25th Workshop on Surveillance and Screening of Blood Borne Pathogens, Tuesday, May 15th 2018, at the Royal Olympic Hotel in Athens

NSS: National Screening laboratory of Sanquin

- Sanquin is the Dutch organization responsible for the blood supply in the Netherlands
- ~340,000 donors (voluntary, non-remunerated donors)
- Annually ~720,000 donations (60% whole blood, 40% plasmaphereses)
- Since 2008 blood screening is centralized at NSS in Amsterdam
 - Blood grouping (AB0, rhD and many others)
 - Serology (HBV, HCV, HIV, HTLV, syphilis)
 - NAT (HBV, HCV, HIV, HEV)

NAT reactives in 2017 (absolute numbers)

- HBV DNA 9
 - HCV RNA 1
 - HIV RNA 3
 - HEV RNA 104 (screening started July 3, 2017)
-
- NAT only's for HBV, HCV and HIV 0

Nucleic Acid Testing laboratory



Cobas 6800 kit controls

- cobas MPX test

Multi target control

- synthetic (armored) HIV-1 Group M RNA encapsulated in MS2 bacteriophage coat protein
- synthetic (armored) HCV RNA encapsulated in MS2 bacteriophage coat protein
- synthetic (plasmid) HBV DNA encapsulated in Lambda bacteriophage coat protein

HIV-1 Group O

- synthetic (armored) HIV-1 Group O RNA encapsulated in MS2 bacteriophage coat protein

HIV-2

- synthetic (armored) HIV-2 RNA encapsulated in MS2 bacteriophage coat protein

Negative control

- negative control: normal human plasma

- cobas HEV test

HEV

- synthetic (armored) HEV RNA encapsulated in MS2 bacteriophage coat protein

Negative control


- negative control: normal human plasma

External quality controls (run controls)


- Specifications
 - Should contain the real viruses (live or inactivated)
 - Have relative low concentrations (challenging for the system)
 - Should be stable when stored below -20°C or -70°C
 - Minimal variation between different lot numbers
 - Calibrated in IU/mL
- To monitor performance of the cobas 6800 robot
 - Over time
 - Detect mechanical issues in an early stage
- To monitor performance of the assay
 - Over time
 - different batches of reagents
- To sleep well
 - Especially when the number numbers of reactive samples is low

The use of run controls at Sanquin

HBV externe controle 2016




HCV externe controle 2016




ThermoFisher MPX controls


HIV-1 externe controle 2016




HIV-2 externe controle 2016




5 IU/mL HBV external control 2017




20 IU/mL HCV external control 2017



75 IU/mL HIV-1 external control 2017



12.5 cp/mL HIV-2 external control 2017



Evaluation of two run controls (BioQControl)


ViraQ Multi-Marker Check 75

- mixture of inactivated HBV, HCV and HIV-1 plasma standards in plasma
- contains 75 copies /mL HBV DNA (gt A), HCV RNA (gt 3a) and HIV RNA (subtype B)
- calibrated against secondary WHO standards: 14 IU/mL HBV, 27 IU/mL HCV, 129 IU/mL HIV
- approx. 5-10 x LOD cobas 6800 MPX test

ViraQ HEV Check 125


- plasma from HEV RNA (gt 3) positive blood donation
- calibrated against secondary WHO standards: 100 IU/mL
- approx. 5x LOD cobas 6800 HEV test

2017: Validation period testing ViraQ Multi-Marker Check 75




HBV run control 2018

14 IU/mL




HCV run control 2018

27 IU/mL




HIV-1 run control 2018

129 IU/mL






HEV run control 2018

100 IU/mL



Comparable Dutch and Belgian[#] data on two reagent lots of cobas MPX on ViraQ Multi-marker Check 75 Control



marker	MPX lot	n NL	n BE	Avg (Std) NL	Avg (Std) BE
HBV	1	9	8	34.28 (0.51)	34.27 (0.42)
	2	42	22	34.33 (0.42)	34.33 (0.31)
HCV	1	9	8	37.71 (0.28)	37.29 (0.55)
	2	42	22	37.43 (0.55)	37.43 (0.64)
HIV	1	9	8	34.32 (0.33)	34.38 (0.34)
	2	42	22	34.21 (0.40)	34.33 (0.31)

[#] data Flemish Red Cross Blood Center, Mechelen, Belgium

Conclusions

- Since 2010 Sanquin has used run controls for the S201 system and cobas 6800 system
 - always very stable results indicating stable test systems and stable run controls
 - problem with a new lot of external controls from ThermoFisher in March 2017
- After 4 months of testing on a daily basis, ViraQ Multi-Marker check 75 and ViraQ HEV check 125 show stable test results (one lot of external controls and two lots MPX/HEV reagents) on the cobas 6800 system
- Comparable Dutch (NSS) and Belgian (Flemish Red Cross Blood Center) results on two reagent lots of cobas MPX on ViraQ Multi-marker Check 75 Control

Acknowledgements



Colleagues from the National Screening Lab Sanquin & Virus Diagnostic Services and colleagues from the Flemish Red Cross Blood Center, Mechelen, Belgium

Analytical sensitivity of Roche PCR assay versions on HBV-DNA standards

HBV genotype A standard	panel	NAT method	n	50% LOD (CI) cp/mL	95% LOD (CI) cp/mL
S0011 VQC-Sanquin WHO 97/746#	PeliCheck	Multiprep Ampliscr.	24	5.9 (3.9-9.3)	72.4 (37.3-211.4)
	Sanquin	Multiprep Ampliscr.	24	3.7 (2.7-5.3)	27.7 (17.6-53.3)
S0011 VQC-Sanquin WHO 97/746#	PeliCheck	Ampliprep Ampliscr.	24	8.6 (5.9-12.9)	88.2 (48.7-215.7)
	Sanquin	Ampliprep Ampliscr.	24	5.3 (3.7-6.9)	35.7 (22.9-69.3)
S0011 VQC-Sanquin WHO 97/746#	PeliCheck	TaqScreen 1.0	12	2.9 (1.6-4.9)	26.9 (13.9-71.0)
	EFS	TaqScreen 1.0	24	3.7 (2.8-5.0)	22.0 (15.5-34.6)
S0010 Eurohep	P0001	TaqScreen 1.0	12	2.3 (1.3-3.8)	14.1 (7.2-56.6)
S0043 BioQ inact.	P0251	TaqScreen 2.0	12	2.8 (1.5-4.3)	23.8 (12.4-99.3)
S0043 BioQ inact. WHO 97/750#	P0031	cobas MPX	12	2.4 (1.4-4.2)	18.6 (9.1-75.9)
	P0023	cobas MPX	12	1.8 (0.93-2.8)	8.0 (4.4-37.4)
S0011 VQC-Sanquin	P0007	cobas MPX	24	1.9 (1.3-2.7)	13.0 (7.7-29.6)

1 IU = 5.33 copies

Analytical sensitivity of Roche PCR assay versions on HCV-RNA standards

HCV genotype 1 standard	panel	NAT method	n	50% LOD (CI) cp/mL	95% LOD (CI) cp/mL
WHO 96/790#	Roche	Multiprep-Ampliscr.	24	19.9 (6.0-27.3)	79.2 (57.3-234)
S0009 VQC-Sanquin	PeliCheck	Ampliprep Ampliscr.	24	16.5 (11.5-22.3)	90.6 (60.4-168)
WHO 96/790#	Sanquin	Ampliprep Ampliscr.	24	7.9 (4.6-15.0)	62.8 (30.0-289)
S0009 VQC-Sanquin	PeliCheck	TaqScreen 1.0	12	4.9 (3.2-7.5)	19.7 (12.0-43.2)
WHO 96/798#	EFS	TaqScreen 1.0	24	6.1 (4.9-8.2)	48.9 (13.4-83.8)
S0109 BioQ gt 3a inact.	P0020	TaqScreen 2.0	12	5.2 (3.3-7.8)	35.2 (19.3-114)
S0109 BioQ gt 3a inact.	P0020	cobas MPX	10	2.5 (1.3-4.4)	15.6 (7.6-77.7)
S0009 VQC-Sanquin	P0019	cobas MPX	60	2.9 (2.3-3.6)	20.4 (14.3-33.3)

1 IU = 2.73 copies

Analytical sensitivity of Roche PCR assay versions on HIV-1 RNA standards

HIV-1 subtype B standard	panel	NAT method	n	50% LOD (CI) cp/mL	95% LOD (CI) cp/mL
WHO 97/656#	Roche	Multiprep-Ampliscr.	24	8.6 (5.1-11.3)	30.4 (23.4-53.4)
S0012 VQC-Sanquin	PeliCheck	Ampliprep Ampliscr.	24	11.8 (8.4-17.2)	76.9 (45.4-182)
WHO 97/656#	Sanquin	Ampliprep Ampliscr.	24	7.0 (5.1-10.1)	72.9 (42.1-164.6)
S0012 VQC-Sanquin	PeliCheck	TaqScreen 1.0	12	4.2 (2.6-6.7)	22.8 (13.0-52.9)
WHO 97/650^	EFS	TaqScreen 1.0	24	4.5 (3.4-5.8)	27.6 (19.1-44.1)
S0041 BioQ inact.	P0251	TaqScreen 2.0	12	2.0 (1.3-2.8)	7.6 (4.9-21.4)
S0041 BioQ inact.	P0026	cobas MPX	12	1.0 (0.6-1.60)	5.8 (3.0-23.2)
WHO 97/650^	P0022	cobas MPX	12	2.7 (1.7-3.9)	5.8 (3.9-24.9)

1 IU = 0.39 copies

^ 1 IU = 0.58 copies

Analytical sensitivity of MPX and HEV test on cobas 6800

Target	95% LOD (IU/mL)	95% CI
HBV	1.4	1.2-1.7
HCV	7.0	5.9-8.6
HIV-1	25.7	21.1-32.8
HIV-2	4.0	3.3-5.2
HEV	18.6	15.9-22.6